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"Truly impressive work! The PROCOR-50's combination of pulsatile flow, energy efficiency, and compact design could significantly expand the applicability of TAHs to small-BSA patients. Innovations like these show how engineering and physiology can merge to create solutions that were previously unimaginable."

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11th BIOMATSEN 2026
International Congress on Biomaterials and Biosensors
Liberty Hotels Lykia / Oludeniz – Turkey
April 2-6, 2026

Characterization of PROCOR-50, a new 50cc pneumatic total artificial heart with low power consumption

Charles-Henri DAVID¹, Said CHABANE¹, Pierre FLEURANCE¹, Emeline MALLET¹, Samuel PLUMEJALUT¹, Thomas SENAGE², Sylvain THUAUDET¹, Jean-Christian ROUSSEL²

1. Procope Medicals, 16 Bd Joliot Curie, 44200 Nantes, France
2. Nantes Université, CHU Nantes, CNRS, INSERM, Institut du thorax, F-44000 Nantes, France


BACKGROUND

- Current total artificial hearts (TAHs) are often limited in patients with small body surface area (BSA), particularly female recipients, highlighting the need for next-generation TAHs suitable for a broader patient population [1,2].
- A new fully implantable TAH – PROCOR-50 – compact and adaptable to different anatomical profiles (Figure 1), is under development and has been pre-validated on cadaveric models (Figure 2). It provides orthotopic biventricular support with two 50 ml ventricles, four mechanical valves, and a closed-loop pneumatic system, including a miniature pneumatic compressor and sequencer technology (Figure 3).
- The objective of this fully orthotopic TAH with intrathoracic technology is to improve quality of life and reduce infection risk from percutaneous cables, while delivering physiological pulsatile flow with low power consumption and full compatibility with transcatheter energy transfer (TET), offering a promising solution for smaller patients.


MATERIALS & METHODS

- Hemodynamic performance of the new TAH was evaluated using a four-chamber mock circulatory loop (Figure 4) reproducing systemic and pulmonary circulations with adjustable resistances and fixed compliances (1 ml/mmHg for the systemic arterial chamber, 10 ml/mmHg for the systemic venous chamber, 2 ml/mmHg for the pulmonary arterial chamber and 5 ml/mmHg for the pulmonary venous chamber) [3]. A water-glycerine test fluid (3.5 cP, 1090 kg/m³) was used to simulate blood properties [4]. Mock pressures (AoP, PAP, LAP, RAP) and TAH ventricular pressures were monitored using pressure transducers, cardiac output was measured via ultrasonic sensors and power with a torque meter and a wattmeter. Device control and parameter tuning were performed through dedicated software.
- Normotensive baseline conditions (RAP = 6 mmHg, PAP = 20 mmHg, LAP = 8 mmHg, AoP = 85 mmHg) were established, rest and exercises were simulated by varying beat rate (90–150 bpm) and cardiac output. Afterload and preload variations were applied to assess ventricular response and overall TAH behavior [5].
- To optimize the hemodynamic profile, preliminary tests were conducted at the laboratories of Procope Medicals [6, 7], leading to the selection of ON-X valves for their optimal balance of Effective Orifice Area (EOA), pressure gradient, and regurgitant fraction (RF): **ON-X valves sizes 25 in mitral and tricuspid positions** (in mitral position at 5 l/min and normotensive conditions: EOA = 2.55 cm², pressure gradient = 2.5 mmHg, RF = 9.94) and **23 in aortic positions** (in aortic position at 5 l/min and normotensive conditions: EOA = 2.86 cm², pressure gradient = 6.70 mmHg, RF = 4.47%).

RESULTS

- The device operates nominally up to a 40 ml stroke volume and achieves a maximum cardiac output of 6 l/min at 150 bpm, with pressures within the normotensive range. The results presented in Figures 5 and 6 correspond to a test at 120 bpm, within the nominal operating range (with 4.8 l/min and a stroke volume of 40 ml).
- Overall power consumption remained below 22 W across the operating range (90 to 150 bpm, 3.5 to 6 l/min) (Figure 7).
- TAH performance is afterload-insensitive, with stroke volume and cardiac output stable despite increased left ventricular pressures and AoP (85 to 135 mmHg) (Figure 8).
- In contrast, with increased preload (leading to higher LAP), stroke volume and cardiac output increase, indicating significant preload dependence (Figure 9).

REFERENCES

[1] HYBRIDHEART Consortium. White paper on requirements and constraints of combining technologies; 2022.
[2] Araba FA. Total Artificial Heart: what metrics should be met for standard care. In: Proceedings of Transcatheter Summit, 2018 Mar 5-8; Phoenix, AZ.
[3] Donnan E. Design of a hydraulic analog of the circulatory system for evaluating artificial hearts. Artif Organs. 1972;6(3):339-349.
[4] Nader F, Skinner S, Brennan M, Kim R, Lammone N, Guller N, Gauthier A, Andrade-Jaramila S, Renaud C, Hardy-Dumontet M-D, Szafrany E, Joly P, Bernhart F, Gervais F. Blood Rheology, Viscosity, Flow, and Flow Rate in Blood Circulation: Rate in Blood Circulation and Effects of Exercise. Frontiers in Physiology. 2019;9:1219. doi:10.3389/fphys.2019.01219.
[5] Cronje M, Dennek M, Truij P, Smit M, Lansen SE, Pothuis ZJ, Buijkoff D, Stehouwer M. Physiological characterization of the SynCardia Total Artificial Heart in a mock circulation system. ASAIO J. 2015;61(11):274-281. doi:10.1089/ASAI.0000000000000102.
[6] ISO 5840-2:2021. Cardiovascular implants—Cardiac Valve Prostheses—Part 2: General Requirements. International Organization for Standardization, Geneva, Switzerland; 2021.
[7] ISO 5840-2:2021. Cardiovascular implants—Cardiac Valve Prostheses—Part 2: Surgically implanted heart Valve Substitutes. International Organization for Standardization, Geneva, Switzerland; 2021.
[8] Guller N, et al. Pothuis M, Kim R. Balancing the ventricular outputs of pulsatile total artificial hearts. Artif Organs. 2023;47(12):1809-1817. doi:10.1111/aor.14641.

CONCLUSION

- The new TAH provides physiological cardiac output up to 6 l/min, with stable hemodynamics and low energy consumption (< 22 W), compatible with TET. The device will be able to balance its left/right flow [8].
- Its compact design makes the PROCOR-50 suitable for small BSA-patients while maintaining physiological pulsatile flow and high efficiency.
- Future work will focus on design optimization, preclinical validation, and the roadmap toward first-in-human implantation.

 PROCOPE
MEDICALS

Merci

PROCOPE-MEDICALS.COM